

# LEWA homogenizer

test system.



# LEWA homogenizer test system. The metering system with an integrated homogenization valve for emulsions.



In an emulsion, droplets that are usually of different sizes are finely dispersed in a second liquid with which they cannot mix. For process engineering reasons, it is often important to reduce the droplet size and make the size distribution as homogeneous as possible. This is necessary in order to make the mixture more stable or to increase the active surface area of the disperse phase. In particular, the LEWA test system can be used to test the aseptic operation of a homogenizer based on a hermetic process diaphragm pump. A downstream spray drying step could thus be supplied directly in a hygienic manner.

### Performance overview:

Discharge pressure up to 400 bar Media flow rate (max. 1,000 cp) up to approx. 150 l/h

#### **Functional principle:**

The flow rate and pressure required to homogenize the process fluid at the two-stage, dynamically controllable homogenization valve are realized by a LEWA triplex process diaphragm pump. After reaching a specified upstream pressure on the suction side via a controlled twin screw pump, the high-pressure diaphragm pump, which is also speed-controlled, first conveys the required flow rate to a preselected pressure P2 of the 2nd stage (valve remains pneumatically closed during this process). Once the pressure and pump capacity have been reached, the homogenization valve of the 1st stage closes until the default value is also reached there. The valve of the 2nd stage is no longer controlled. The ratio is usually P2 <  $1/5 \times P1$ , but can be freely selected. The now leading homogenization pressure of the 1st stage is also maintained by pneumatic valve control of the dispersion gap. The system can be operated both standalone (batch operation) in a circuit and continuously integrated into the system. Fluid is supplied either via the heatable vessel on the system side or via a separate media connection.


#### Essential components of the LEWA test system:

- LEWA triplex process diaphragm pump with diaphragm rupture indicator, incl. motor 5.50 kW
- Twin screw pump with double-acting purged mechanical seal, incl. motor 2.2 kW
- Equipment cabinet (at the rear of the system), hygienic design
- Heatable double-walled vessel, 60 l, unpressurized, lockable, max. 110 °C
- Two-stage sterilizable homogenization valve, pneumatically controlled, max. 400 bar/150 l/h
- Diaphragm valves, flows switchable depending on the operating mode
- Safety valve (automatically cleanable inline) to protect the system against overpressure
- Various pressure and temperature transmitters
- Warning elements: Signal lights and horn
- Water connection with piping for cooling the control cabinet and flushing the mechanical seal of the twin screw pump

The above-mentioned components are piped hygienically and mounted on a common base frame. The system is suitable for washdowns. The specially developed LEWA smart control controller is integrated into a control cabinet mounted on the metering frame. The "Automatic" and "Manual" operating modes can be selected in the control system. Automated processes such as CIP, SIP or emptying are stored in step chains.



# Flow diagram of the LEWA homogenizer test system



# Dimensions and connection data

# Dimensions:

- 2.56 x 1.66 x 2.16 m (L x W x H) plus swing area of the doors
- Weight: approx. 1,800 kg, filled

## Electrical connection data:

- Voltage: 400 V/50 Hz
- Performance: 12.5 kW

# Other connection data:

Compressed air and cooling water:
For details see below (connection list)

# **Connection list**

Number Size	Pressure rating	Туре	Description for operating mode	
			Standalone mode	Within a process line
DN25	10 bar	DIN 11864-1 threaded connection	Drain suction line	Product inlet
DN25	10 bar	DIN 11864-1 threaded connection	Drain/cleaning pressure line	Product outlet
Ø 12 mm		Pipe end	Cooling water outlet control cabinet/ screw spindle pump	
DN25	1 - 5 bar	DIN 11851 threaded connection	Cooling water inlet control cabinet/ screw spindle pump	
Ø6mm			System-internal pneumatic line	
Ø6mm			System-internal pneumatic line	
Ø8mm	2 - 3 bar	Compression fitting	Inlet cleaning (CIP)/sterilization (SIP) Homogenization valve	
Ø8mm	2 - 3 bar	Compression fitting	Outlet cleaning (CIP)/sterilization (SIP)	
	6 - 10 bar	Quick release coupling	Compressed air supply (quality class 4.4.4 in accordance with ISO8573-1)	
	DN25 DN25 Ø 12 mm DN25 Ø 6 mm Ø 6 mm Ø 8 mm	DN25     10 bar       DN25     10 bar       Ø 12 mm	DN2510 barDIN 11864-1 threaded connectionDN2510 barDIN 11864-1 threaded connectionØ 12 mmPipe endDN251 - 5 barDIN 11851 threaded connectionØ 6 mmImage: Compression fittingØ 8 mm2 - 3 barØ 8 mm2 - 3 barCompression fitting	SizePressure ratingTypeStandalone modeDN2510 barDIN 11864-1 threaded connectionDrain suction lineDN2510 barDIN 11864-1 threaded connectionDrain/cleaning pressure lineØ 12 mmPipe endCooling water outlet of screw spindle pumpDN251 - 5 barDIN 11851 threaded connectionCooling water inlet co screw spindle pumpØ 6 mmDin 11851 threaded connectionSystem-internal pneumØ 6 mmCompression fittingInlet cleaning (CIP)/ste Homogenization valveØ 8 mm2 - 3 barCompression fittingOutlet cleaning (CIP)/ste Homogenization valveØ 8 mm2 - 3 barCompression fittingOutlet cleaning (CIP)/ste Homogenization valve



Your local representative:

**LEWA GmbH** Ulmer Str. 10 71229 Leonberg